

DESCRIPTION
INSERTER, IMPLANT
BONE GRAFT PACKING BLOCK
BONE GRAFT TAMP
TAMP, IMPLANT
END PLATE RASP
QUICK RELEASE HANDLE
TRIAL, 14 X 12 mm X (5-12 1mm incr.) 0°
TRIAL, 15 X 13 mm X (5-12 1mm incr.) 0°
TRIAL, 14 X 12 mm X (5-12 1mm incr.) 6°
TRIAL, 15 X 13 mm X (5-12 1mm incr.) 6°



Phone: 760-295-3600  
 Fax: 760-295-3610  
 Address: 3270 Corporate View, Suite A  
 Vista CA 92081  
 Email: [customerservice@nextorthosurgical.com](mailto:customerservice@nextorthosurgical.com)



**EZ-TRACK™**  
 CERVICAL INTERBODY CAGE SYSTEM



**Axial Window:**  
Axial window for the acceptance of autogenous bone graft.



**Tracking Ridge:**  
Provides precise insertion direction, advanced control and safety.



**Biconvex:** Contoured to match the vertebral endplates.



FOOTPRINTS	WIDTH (mm)	DEPTH (mm)	HEIGHT (mm)	LORDOSIS TYPE
12 X 14	14	12	5 - 12 (1mm incr.)	0°,6° FLAT, DOMED
13 X 15	15	13	5 - 12 (1mm incr.)	0°,6° FLAT, DOMED
14 X 16	16	14	5 - 12 (1mm incr.)	0°,6° FLAT, DOMED

**FLAT**

**DOMED**

