


DESCRIPTION
AXIAL HANDLE, LARGE
T-HANDLE, LARGE
POSTERIOR PADDLE SHAVER, SIZES 5-14mm
PLIF, TRIAL, 9 x 24 x (7-14mm) -5°
OTLIF, TRIAL, 9 x 24 x (7-14mm) -0°
TLIF, TRIAL, 11 x 28 x (7-14mm) -5°
POSTERIOR PADDLE DISTRACTOR, SIZES: 7-14mm
PLIF, TAMP, STRAIGHT, IMPLANT
PLIF, TAMP, ANGLED, IMPLANT
PLIF, BONE PACKING BLOCK
PLIF, BONE GRAFT TAMP
PLIF, INSERTER, STRAIGHT
PLIF, INSERTER, PISTOL
PLIF, MALLET
PLIF, CANCELLOUS BONE FUNNEL, TUB
PLIF, CANCELLOUS BONE FUNNEL, SHAFT
PLIF, CANCELLOUS BONE FUNNEL, PUSHER
PLIF, SLAP HAMMER

Phone: 760-295-3600
 Fax: 760-295-3610
 Address: 3270 Corporate View, Suite A
 Vista CA 92081
 Email: customerservice@nextorthosurgical.com

EZ-TRACK™
 POSTERIOR LUMBAR INTERBODY SYSTEM



Bone Trap:

The Bone Trap feature retains bone graft throughout insertion.


Tracking Ridge:

Provides precise insertion direction, advanced control and safety.



Biconvex: Contoured to match the vertebral endplates.

CAGE TYPE	WIDTH (mm)	LENGTH (mm)	HEIGHT (mm)	LORDOSIS
PLIF	9, 11, 13	20 - 30 (2mm incr.)	7 - 14 (1mm incr.)	5°
TLIF (5°)	9, 11, 13	22 - 32 (2mm incr.)	7 - 14 (1mm incr.)	5°
TLIF (0°)	11	28	7 - 14 (1mm incr.)	0°
OTLIF	9, 11, 13	20 - 34 (2mm incr.)	7 - 14 (1mm incr.)	0°



PLIF

TLIF

OTLIF